

# BIOTHESIA COUPLED IMAGING SYSTEM

## Specifications

Make & Model	VILBER FUSION FX system
Analysis	Chemiluminescence, Bioluminescence, UV fluorescence, visible imaging.
Imaging	DARQ-7 camera- 4.2 megapixel , HSR technology.
Sensitivity	V.084 lens for very faint light condition
Filter	590 nm emission maxima from gels stained with EtBr, Gel Red, Gel Green, Lumitein, SYPRO Ruby
Applications	Western blot, small animals, DNA gel, Protein gel imaging.