

ISOELECTRIC FOCUSING SYSTEM (IEF)

Specifications

Make & Model	Bio-Rad PROTEAN [®] i12 [™] IEF System Focusing tray
Material	Polycarbonate
IPG strip range	pH 3-10 ReadyStrip [™]
IPG strip length (cm)	7, 11 & 17
Channel volume	7 cm – 7 ml, 11 cm – 10 ml 17 cm – 14.2 ml
Capacity	1-12 IPG strips / tray Rehydration/equilibration tray
Material	Polystyrene
IPG strip range	pH 3-10 ReadyStrip [™]
IPG strip length (cm)	7,11 & 17
Channel volume	7 cm – 6.8 ml, 11 cm – 96 ml 17 cm – 14.2 ml
Capacity	1-12 IPG strips / tray
Applications	Protein profiling, Biomarker discovery, Functional proteomics.